Know What's in Your Membranes: Quantification of Factors in Placental Allografts³

A variety of membrane products have been developed for applications focus on exploiting the available growth factors and adding strength to AM and AC. These include lamination, cryopreservation method results in the optimization of different properties of the membrane. Although there is a vast array of information collected on amnion/chorion membranes, the studies often utilize enzymatically digested, macerated and/or centrifuged extracts. While using the digested, pre-processed membrane to release all proteins may give us a clear look at the potential of the membrane, it does not allow us to determine the factors available to a patient. For treatment with membrane product, the tissue is usually sterile, intact and laid on a wound or treatment area. The factors present as available to the treatment area should be quantified and reported in a standardized manner. To address this, we characterized AM and AC from the perspective of factors eluted from a specific area (per square centimeter) of finished product, as would be eluted to patient tissue. Analysis of dehydrated, sterilized amnion/chorion revealed elution of factors that are conducive to wound healing.

Methods **A Standardized Process**

- Donated human placentas were acquired and processed in accordance with FDA Good Tissue Practices and AATB standards.
- AM or AC was isolated from the placenta and BioREtain[®] processed, resulting in sterile AMⁱ or ACⁱⁱ.
- The final sterilized products were used for all tests.
- 8 mm biopsy punches were taken from tissue.
- Punches were placed in 500 uL DPBS^{Ca-Mg-} at 37° C for 72 hours, rocking.
- The supernatant was collected, clarified by centrifugation and stored at -80° C until use.
- Custom 5-plex cytokine panels and ELISA assays were performed.
- Results were plotted using a 5-parameter logistic standard curve.
- The average pg/ml of growth factor was multiplied by the milliliters of eluate= total GF content in the eluate.
- The total GF content in the eluate was divided by 0.503 (8mm punches have an area of 0.503 cm²) to reveal the pg of growth factor in 1 cm^2 of tissue.

 $\frac{GF\left(\frac{pg}{ml}\right)*ml}{0.503 \ sqcm} = pg \ of \ growth \ factor \ per \ cm^2.$

This can be utilized for any tissue, volume of eluate used and any sample punch size.

^{i.}Vendaje[®] ^{ii.} Vendaje AC[®]





Taylor J Sabol, MS^{1,2}, Jason V Matuszewski, BS¹, Wendy W Weston, PhD^{1,2}

¹ BioStem Technologies, Inc., ² CTBS

Abstract

